

**B. Graczyk, Inc.**  
**Universal™ / Surgical Eye Care Products**

Product Overview Sheet

Name: Universal™ Concentrated Surgical Instrument Cleaner and Lubricant

1. FDA Approved
2. Laboratory Tested
  - Antimicrobial / Rinses Residue Free
  - No Cytotoxicity / No Cell Lysis, Intracytoplasmic Granules
  - Synthetic and Mineral Formulation / No Animal Bi-Products
3. Clinically Proven to be Safe and Effective
4. Approved by all major microkeratome and surgical instrument companies

Part No. Description:

40-803 120 ml Calibrated Plastic Bottle / Cleans Approximately 100 LASIK Instrument Sets  
40-805 240 ml Calibrated Plastic Bottle / Cleans Approximately 200 LASIK Instrument Sets  
40-810 960 ml Calibrated Plastic Bottle / Cleans Approximately 800 LASIK Instrument Sets

Featured Benefits:

1. Multiple Cleaning Applications - Washer Sterilizers, Ultrasonic Cleaners, Manual / Hand Washing and Soaking Instruments.
2. Concentrated Cleaner, 30 ml to one 3.8 Liter of Distilled or Deionized Water.
  - Universal™ Concentrated Cleaner and Diluted Universal™ Cleaner have a stable shelf life for up to 5 years.
3. Economical to use - as low as 0.4 cents per cleaning procedure.
4. Eliminates the need to stock additional cleaning and lubricating products.

Special Benefits:

1. Citrus base aids in removing residues and biofilms.
2. Emulsifiers to soften dried blood, protein deposits and other surface contaminants.
3. Conditioners for hard water.
4. Stabilizers to help retard scale buildup.
5. Advanced Penetrating Lubricant that bonds to the surfaces of the metal.
6. Rinses Residue Free.

FEATURES:	ADVANTAGES:	BENEFITS:
Combination Cleaner & Lubricant	Eliminates dual processing	Cost Saving - Time Saving
Mixes / Dilutes Instantly	Will not clump	Leaves No Residues
Lubricant Bonds to surfaces of metals	Deep penetration without buildup	Rinses Residue Free
Removes blood, protein and stains	Minimal hand washing is required	Instruments are clean and work freely

Application:

Provides a safe, effective and economical method of meticulously cleaning surgical instruments.

Use:

Dilute 30 ml of the Universal™ Cleaner into one 3.8 Liter of distilled or deionized water; use the diluted cleaner until the container is empty. Universal™ Cleaner has a stable shelf life for up to 5 years.

Projected Costs Per Procedure:

- 40-810 Universal™ Concentrated Surgical Instrument Cleaner and Lubricant **960 ml Bottle**
- 30 ml of Concentrated Cleaner to one (1) 3.8 Liter Demineralized / Filtered water.
  - 960 ml Concentrated Cleaning Solution Makes 32 3.8 Liters of Diluted Cleaner.
  - Cost Per 3.8 Liter in Use = \$ 1.50
  - Cleaning Procedures Per Diluted Liter in Use = 25 Instrument Sets
  - Cleaning Procedures for 32 Diluted Liters in Use = 800 Instrument Sets
  - Cost Per Cleaning Procedure = \$ 0.04 / Economical to Use